

Dental Equipment

Retail Price List

April 2017



Standard Warranty Policy

All Belmont products are guaranteed against mechanical, electrical and workmanship faults, under normal working conditions, for a period of one year from date of delivery. Paintwork, plating, natural wear/tear of all tubings are not warranted unless they are proved to be defective at time of delivery. OEM factory built-in options (i.e. scalers, motors and CVS kits) carry a warranty as provided by the individual manufacturer.

FREE Extended Warranties

On condition that the purchaser returns the *Extended Warranty* form to Takara Belmont (UK) Ltd, the above *Standard Warranty* will be extended to ~

5 Years

for all Belmont Dental Chairs, Dental Operating Lights, Cleo Dental Units, tbCompass Dental Units and Voyager Dental Units.

2 Years

for Intra-Oral X-Rays.

(A detailed warranty policy is available upon request.)

PLEASE NOTE: ISO 13485 Certification

All of our dental products are supplied with a full CE mark and certification to the standard complying with ISO 13485. Any modifications or alterations made after certification by us may affect your warranty as the responsibility for CE and electrical safety will transfer to the modifier.

Intra-Oral X-Ray Units

097 BELRAY II (universal compatibility)

097WK Belray II

wall-mounted, with 80cm primary arm, wall plate, scissor-type arm

£2,805

PHOT-X IIS DC X-RAY (universal compatibility)

505WK

wall-mounted, with 80cm primary arm, wall plate, scissor-type arm

£3,290

PHOT-X II DC X-RAY (universal compatibility)

303FM

mobile base

£3,615

303CK

ceiling-mounted, with ceiling-mount pole & bracket

£3,980

303RK

room type, with stool

£4,090

303FK

floor-mounted base

£3,080

097 Belray II

features include:

- 70 KVP
- 0.7mm focal spot
- Digital time display
- Digital imaging compatible
(min. exposure 0.02 sec)
- Manual override for total flexibility
- Touch control for 2 film settings & digital imaging
- Lockable touch pad control
- Optional cone length & tube current
- Rectangular Collimator
- Hand exposure switch assembly

models include:

097WK Wall-mounted, with 80cm primary arm, wall plate, scissor-type arm (max. reach 185cm)

Phot-X IIS DC X-Ray

features include:

- 60 or 70 KVP options
- 0.4mm focal spot
- Digital time display
- Radiation dose display
- Compatible for digital imaging
(min. exposure 0.01 sec)
- Manual override for total flexibility
- Touch control for 2 film settings & digital imaging
- Lockable touch pad control

- Optional cone length, tube current & tube voltage
- Rectangular Collimator
- Hand exposure switch assembly
- Error code history recorded

models include:

505WK Wall-mounted, with 80cm primary arm, wall plate, scissor-type arm (max. reach 192cm)

Phot-X II DC X-Ray

features include:

- 60 or 70 KVP options
- 0.7mm focal spot
- Digital time display
- Compatible for digital imaging
(min. exposure 0.01 sec)
- Manual override for total flexibility
- Touch control for 2 film settings & digital imaging
- Lockable touch pad control
- Optional cone length, tube current & tube voltage
- Rectangular Collimator
- Hand exposure switch assembly

models include:

- 303FM** Mobile base
- 303CK** Ceiling-mounted, with ceiling-mount pole & bracket
- 303RK** Room type, with stool
- 303FK** Floor-mounted base

Optional extras for 097 Belray II

£FOC	30cm primary arm in lieu of 80cm
£65	100cm primary arm in lieu of 80cm
£260	Long cone for paralleling technique

Optional extras for Phot-X IIS

£FOC	30cm primary arm in lieu of 80cm
£65	100cm primary arm in lieu of 80cm
£260	Long cone for paralleling technique
£265	Operating Light bracket for 303CK

Standard Colours

Smooth-to-touch for easiest hygiene maintenance



For Seamless upholstery



Black
CG5



Light Grey
CG3



Beige
CG8



Ivory
CG2



Light Green
CG23



Teal
CG18



Blue
CG19



Regimental Blue
MR1



Crocus
MR2



Mauve
CG16



Pink
CG22

See product images in your preferred colours
and request swatches online at belmontdental.co.uk

IMPORTANT NOTE:

Cancelled or amended orders of Ultrasoft upholstery will incur the following fee:
£150 per Luxury chair, £75 per Seamless chair, and £50 per stool.

Ultrasoft Colours

Soft-to-touch for superior comfort



For Luxury and Seamless upholstery



Lagoon
291-2527



White
291-0001



Ivory
291-3700



Sand
291-3850



Pecan
291-3778



Taupe
291-3779



Sassy
291-5250



Fudge
291-3666



Apricot
291-8243



Wet Sand
291-3491



Hide
291-3610



Salmon
291-6588



Hot Sauce
291-8241



Poppy
291-1347



Toro
291-1334



Red
291-1176



Sorbet
291-6580



Chianti
291-1312



Hot Chocolate
291-3205



Pine Cone
291-3183



Sonora
291-4454



Orchard
291-4262



Parrot
291-4460



Eucalyptus
291-4430




291-5243
Lemon



Cyan
291-2553



Bayou
291-2685



Silver Pearl
291-5668



Riviera
291-2552



Dove Grey
291-5666



Baltic
291-2551



Admiral
291-2556



Plum
291-9385




Stone
291-3602



Schooner
291-5806



Granite
291-5708



Raven Wing
291-5607

 **Belmont**

Takara Belmont (UK) Ltd.

London Head Office & Showroom

Belmont House, One St. Andrews Way, Bow, London E3 3PA

T. +44 (0) 20 7515 0333

F. +44 (0) 20 7987 3596

Manchester Showroom

7 Quays Reach, Carolina Way, Salford, Manchester M50 2ZY

T. +44 (0) 161 743 9992

E. dental@takara.co.uk

All showroom visits by appointment only

All prices listed are shown exclusive of VAT.

The manufacturers reserve the right to alter the prices and specifications without prior notice.

 **Belmont**

belmontdental.co.uk